

10 MELAG Vacuklav 44-B

15 Program: Universal program
20 Program type: 134°C wrapped 5:00 min
25 Date: 01.11.2005
30 Daily batch: 02 Total: 00368

40 Universal program successfully ended
42 =

45 Temperature: 135.3 +0.21/-0.39°C
50 Pressure: 2.16 +0.08/-0.22 bar
55 Plateau time: 05 min 00 s
60 Conductivity: 2 µS/cm (226:18.2)
65 Start time: 13:26:39
70 End time: 13:57:39 (31:00 min)

80 SN:200344-B1001 BO V2.1
>> never change code in following row <<
0D000008005302856BEA030303E0037531583474
>> Proof of veracity batch log files <<

Step	Time t[m:s]	P[mbar]	T[°C]
SK11	0:16	0:16	1603 99.2
SK12	0:39	0:23	1284 104.7
SK11	0:52	0:13	1614 107.5
SK12	1:16	0:24	1286 111.4
SK21	1:21	0:05	1624 111.6
SK22	1:45	0:24	1284 116.4
SK21	1:51	0:06	1616 116.2
SK22	2:15	0:24	1281 118.9
SK21	2:21	0:06	1630 118.4
SK22	2:44	0:23	1290 120.0
SF12	3:04	0:20	478 115.9
SF13	3:32	0:28	1622 117.6
SF21	3:40	0:08	1290 118.5
SF22	4:24	0:44	177 109.6
SF23	5:05	0:41	1814 117.6
SF31	5:17	0:12	1287 118.7
SF32	6:01	0:44	195 110.0
SF33	6:44	0:43	1926 118.6
SF41	6:57	0:13	1290 119.6
SF42	7:27	0:30	382 111.3
SF43	7:50	0:23	1737 116.6
SH01	8:34	0:44	2702 129.7
SH02	8:54	0:20	2837 131.6
SS01	9:41	0:47	3075 134.1
SS02	14:41	5:00	3171 135.4
SA00	15:18	0:37	1293 110.7
ST01	19:18	4:00	60 79.2
ST02	21:18	2:00	661 94.8
ST03	25:19	4:01	59 90.7
ST04	26:20	1:01	662 94.6
ST05	30:20	4:00	64 86.8
SB10	30:38	0:18	801 90.7
SB20	30:59	0:21	1001 92.6
SP-E	30:59	0:00	1001 92.6

MELAdoc
The documentation system



*I want to have **documented** safety,
for our patients and for
our medical practice!*

Quality - Made in Germany

MELAG
Evidence Based Sterilization

The components of documented instruments

1. Process security

The MELAG autoclaves *Vacuklav*[®], *Vacuquick*[®], *Euroklav*[®] and *Cliniklav*[®]25 are equipped with an advanced monitoring system. This system monitors those parameters that are critical for the success of sterilization: i.e., pressure, temperature and time. High-quality precision sensors provide highly exact measured values.

2. Batch monitoring for Class B

MELAcontrol[®] is a test system conforming to EN 867-5 that serves for batch monitoring and for testing the function of the fractionated pre-vacuum of a *Cliniklav*[®]25 or of a "Class B" autoclave. The *MELAcontrol*[®] unit consists of a test receptacle (the so-called helix) and 250 indicator strips. *MELAcontrol*[®] simulates and tests the steam penetration of a long, narrow-bore hollow body.




3. Batch documentation

The log printer *MELAprint*[®]42 records the process parameters for the sterilization cycle. This printer can also be connected to the serial interface of older MELAG models that already have a serial interface.

The *MELAwir*[®] software offers an alternative without paper. It is installed on a PC in the doctor's practice: The PC is directly connected to the autoclave (*Vacuklav*[®]30-B, 31-B, 23-B, 24-B, *Euroklav*[®] and *Cliniklav*[®]25) via the serial interface. Instead of sending data to the *MELAprint*[®]42, this configuration automatically transfers all data to the PC in the practice, where it is stored. In addition, this solution allows entry of additional important information: e.g., load details, the name of the operator, etc.

The new premium autoclaves *Vacuklav*[®]40-B/41-B/43-B/44-B and *Vacuquick*[®]13-B/14-B do not need the *MELAwir*[®]-software. They are intended for direct integration into networks in medical practices, via a built-in Ethernet interface. In addition, it is possible to store the logs on a *MELAflash* CF Card. This makes it possible to keep sterilization records without connecting the autoclave to external devices. A card reader acquires the records from the CF Card and transmits this recorded data to the PC. A standard editor then enables display of the records on the PC.

A network adapter is required to connect the *MELAprint*[®]42 logger to a premium autoclave.



ment preparation

MELAdoc – the system that truly closes a gap

Many doctors' practices document sterilization processes only in a daily log. The MELAdoc Label Printer facilitates and simplifies batch documentation, and enables tracing the data for a certain batch.

This solution makes it easy to trace back and match up a particular instrument or a certain patient to all essential data from the sterilization process: i.e., the date of sterilization, the batch number, the person providing approval clearance for use of the instruments, the sterilizer that was used, and the sterilization expiration date.

a. Documentation of a batch

The MELAdoc documentation system prints a label after each sterilization process, and attaches this label to the Sterilization Journal for documentation of the batch.

b. Designation

After sterilization, MELAdoc affixes a label to designate that the sterilization packaging with

the sterilized instruments has in fact been fully satisfactorily processed. This procedure satisfies all the conditions for enabling the operator of the systems to correctly provide approval clearance of a sterilized batch. The MELAdoc system avoids all the dangers of older systems: e.g., puncture of the packaging by ball-point pens, or partial dissolving of the plastic

packages (with resulting microperforation) by permanent markers

c. For surgical records

MELAdoc labels are multiple-adhesive. This means that they can be peeled off from the instrument packaging after use, and simply attached to the patient's records.

d. Simple tracing

In case that a patient makes a claim for damages for secondary infection, the MELAdoc labels in the patient's records provide all information on the correct sterilization process for the instruments used. The labels allow tracing the patient's data back to the responsible operator and to the equipment used (sterilizer number). The information on the sterilization data and the batch number enable tracing the process back to the sterilization logs (MELAprint®42, MELAflash CF Card, MELAwin® etc.). This procedure provides the required proof: that the instruments were properly sterilized.

Staff number	Sterilizer number	Batch number
02	1	11
Sterilized on:		
07.05.2005		
use before:		
07.11.2005		




Is documentation necessary?

In case of litigation

If a patient institutes legal proceedings for recovery of damages against a medical institution, in the belief that he or she had been harmed (for example, after use of non-sterile instruments), long-standing legal interpretation in many legal systems requires that the patient must prove that the attending physician has committed malpractice and has thereby harmed the patient. As a rule, however, the

patient will be unable to produce the required evidence: i.e., in most cases he or she will not be able to prove his claims or collect for damages. Jurisprudence in many countries, consequently, has institutionalised the practice of reversal of the burden of proof here. This means that the physician's practice or the medical facility must prove that it has not committed malpractice.

If the attending physician can present a complete and correct record of technically satisfactory preparation of the sterilized instruments, established judicial precedent in many countries includes the ruling that the recorded measures have in fact been carried out, and that no liability exists



Quality and precision for the most exacting in hygienic requirements

More than 50 years ago, MELAG began in Berlin to specialise in manufacture of sterilization equipment. Verification of its great, international success has been the sale of more than 355,000 MELAG units. Decades of experience, modern computer-aided production in the Berlin plant, employ-

ment of highest-quality materials, and experienced and reliable staff – all have made MELAG equipment renowned for quality and convenience.

MELAG corporate philosophy includes systematic specialisation in hygiene and sterilization, in a limited line of products.

With the equipment and accessories it produces for legally stipulated equipment preparation, flawless sterilization, and safe sterile storage of sterilized articles, the entire MELAG staff is proud to have made an essential contribution to the protection of the health of patients and of medical staff.



MELAG

www.melag.com

11.05.5E-GD-MDOC