



Bowie & Dick Type Autoclave Test Pack

INTERSTER 



Bowie & Dick Type Autoclave Test Pack



After an intense development process Interster International has developed a new Bowie & Dick Type Autoclave Test Pack which meets all requirements of the ISO 11140-4. The Bowie & Dick Type Autoclave Test Pack is a Type 2 chemical indicator packed in special paper with a cardboard outer packaging which imitates the original Bowie & Dick test. The Bowie & Dick Type Autoclave Test Pack tests if steam penetration in the autoclave is at sufficient level through colour change of the chemical indicator. The change of colour of the chemical indicator will turn from blue to light pink and indicates if all parameters (time, pressure and steam penetration) are sufficient to approve the process.

Material

The Bowie & Dick Type Autoclave Test Pack consists of a chemical indicator which is packed in between two layers of special paper. The chemical indicator and special paper are packed in a cardboard Lean outer packaging which can be opened very easily after usage. The indicator contains a specially developed indicator ink which is affixed through a silkscreen printing process. The product does not contain heavy metals.

Advantages


The new Bowie & Dick Type Autoclave Test Pack has the size of 102 x 76 x 30 mm. Due to the fact the Bowie & Dick Type Autoclave Test Pack has a smaller size, it is possible to pack more efficiently. This has resulted in a very economic transport packaging.

Bowie & Dick Type Autoclave Test Pack has been tested under supervision of the BSI Assurance UK Ltd. The test pack was positively certified with a BSI Kitemark (KM 633144) and completely meets the ISO 11140-4 on all possible Bowie & Dick cycles (B1, B2 and B3).




BS EN-ISO 11140-4:2014
KM 633144

As a result of the starting colour of the Bowie & Dick Type Autoclave Test Pack being blue and the final colour needed being light pink, a clear contrast is visible between a successful and unsuccessful process. This eliminates any doubt between a before and after result of the process. On the image below an unused and used indicator sheet are displayed, the same way they are packed within the Bowie & Dick Type Autoclave Test Pack.



Before use



Correct color change*

*This image **cannot** to be used as reference colour for assessing Bowie & Dick Type Autoclave Test Pack indicator.



Article number	3FSKS610860
Description	Bowie & Dick Test Pack
Unit (box)	40 pieces
Ster. time/temp.	3,5 min. - 134 °C

Interster International B.V.

Samsonweg 2, 1521 RC Wormerveer

T: +31 75 627 22 41 - E: export@interster.nl

www.interster.nl/exporange