

ISP® Sterilization Environment Utility Marker

INTERSTER 




ISP® Sterilization Environment Utility Marker




For marking sterilization packaging material, with special autoclave low odour ink, without additor or oxylene/ toluene. The ink is steam, gamma, formaldehyde, plasma and EO resistant.


The pen is designed for light touch writing to avoid damage to the packaging material. The solid writing point of the marker ensures a constant writing quality and long durability. The ink is dry immediately and does not fade.



Article number	3FSZB710824
Description	Marker
Unit (box)	25 pieces
Colour	Black



Article number	3FSZB710826
Description	Marker
Unit (box)	25 pieces
Colour	Blue



Article number	3FSZB710828
Description	Marker
Unit (box)	25 pieces
Colour	Red